

SYSTEM COMPONENTS:

414-00-300 - NEXA™ DEVICE KIT

- NEXA DEVICE
- NEXA CARRY CASE
- NEXA POWER SUPPLY with region specific adaptor
- INSTRUCTIONS FOR USE

414-00-501 - NEXA™ DRESSING (BOX OF 5)

- NEXA BLACK RETICULATED FOAM 150mm x 75mm x 20mm
- NEXA PRESSURE CONTROL CONNECTOR, DRESSING PAD AND TUBING
- POLYURETHANE DRAPE

414-00-502 - NEXA™ FLUID CONTAINER PACK (BOX OF 5)

- 125ML FLUID CONTAINER & NEXA® DYNAMIC FLUID PUMP AND TUBING

Nexa Medical Limited
Suite 5-6 Willow View, Southfields,
Church Lane, Bournemouth, BH22 8TR.
UK www.nexamedical.com
info@nexamedical.com



414-10-001 (A)

MIN MEDICAL

MIN Medical GmbH
Franz Schalk Platz 9/2
1130 Wien
Austria

Mobil: +43 660 107 67 67
Fax: +43 1 877 68 29
Web: www.minmedical.at
Email: office@minmedical.at



NEXA™

THE NEXADVANTAGE FOR
NEGATIVE PRESSURE WOUND THERAPY





OVERVIEW

The NEXA Negative Pressure Wound Therapy (NPWT) system has been developed to provide traditional, clinically proven therapy for mobile patients with chronic wounds.

With three new innovations in NPWT, NEXA offers more opportunity to transition patients from hospital to the community and reduce the risk of further re-admission due to complications with their wounds.

WHAT IS THE NEXADVANTAGE?

Conventional NPWT devices use a vacuum air pump to generate negative pressure which means exudate removal is dependent on the leak rate of the dressing.

NEXA utilises a Dynamic Fluid Pump which improves exudate removal and provides consistent negative pressure thanks to the novel Pressure Control Connector located close to the wound. The performance of NEXA is also not affected by height differences between the device and the wound. Exudate is drawn into a flexible bag, entirely separated from the pump mechanism.



The NEXA DEVICE

- One button operation for simplicity and ease of use.
- 6 week operational life and approved for multiple patient use.
- 10 hours battery life and 4 hour recharge time.
- Lightweight and small, suited to mobile patients.
- Carry case supplied to enable mobile use of the system.



The NEXA Dynamic Fluid Pump

is constantly active, removing exudate from the wound irrespective of vacuum levels, thereby reducing the potential for skin maceration.



The NEXA FLUID CONTAINER PACK

- 125ml flexible fluid container.
- NEXA Dynamic Fluid Pump.
- Long tube length of 1 metre when connected to the NEXA DRESSING



The NEXA DRESSING

- Sterile, single use dressing set.
- **NEXA Pressure Control Connector**, maintaining constant negative pressure, determined by the dressing selection
- **NEXA Black Reticulated Foam**, attributed to effective negative pressure wound healing, as supported by the clinical evidence.

