

PreludeSNAP™

Splittable Sheath Introducer

The Prelude SNAP from Merit Medical features a **SEAMLESS** sheath-to-dilator transition to reduce vessel trauma and provide ease of insertion through tortuous anatomy. The Prelude SNAP is also **EASIER** to break, requiring 66% less break force than its predecessor, making it effortless to remove and helping ensure leads stay in place.




PRODUCT FEATURES

Robust Valve Design

Tight seal to protect against air embolism and blood loss.

Improved Grip

Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.



Easier to Break

Requires 66% less break force than its predecessor.*

Seamless Transition

Smoother sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.

KIT CONTENTS

- Splittable Sheath
- Dilator
- Guide Wire 0.038" x 50 cm or 0.038" x 80 cm
- Merit Advance Needle 18G x 7cm
- 12 mL Syringe



Competitors' sheath-to-dilator transitions



Our seamless sheath-to-dilator transition

ORDERING INFORMATION

Catalog Number with Sideport	French Size	Sheath Length	Hub Color
PLS-1006	6	13 cm	Green
PLS-1007	7	13 cm	Orange
PLS-1008	8	13 cm	Blue
PLS-1008.5	8.5	13 cm	Light Blue
PLS-1009	9	13 cm	Black
PLS-1009.5	9.5	13 cm	Gray
PLS-1010	10	13 cm	Fuchsia
PLS-1010.5	10.5	13 cm	Light Fuchsia
PLS-1011	11	13 cm	Yellow
PLS-1012	12	13 cm	Brown
PLS-1012.5	12.5	13 cm	Light Brown
PLS-2506	6	25 cm	Green
PLS-2507	7	25 cm	Orange
PLS-2508	8	25 cm	Blue
PLS-2509	9	25 cm	Black
PLS-2510	10	25 cm	Fuchsia
PLS-2510.5	10.5	25 cm	Light Fuchsia
PLS-2511	11	25 cm	Yellow

Catalog Number without Sideport	French Size	Sheath Length	Hub Color
PLSX-1006	6	13 cm	Green
PLSX-1007	7	13 cm	Orange
PLSX-1008	8	13 cm	Blue
PLSX-1008.5	8.5	13 cm	Light Blue
PLSX-1009	9	13 cm	Black
PLSX-1009.5	9.5	13 cm	Gray
PLSX-1010	10	13 cm	Fuchsia
PLSX-1010.5	10.5	13 cm	Light Fuchsia
PLSX-1011	11	13 cm	Yellow
PLSX-1012	12	13 cm	Brown
PLSX-1012.5	12.5	13 cm	Light Brown
PLSX-2506	6	25 cm	Green
PLSX-2507	7	25 cm	Orange
PLSX-2508	8	25 cm	Blue
PLSX-2509	9	25 cm	Black
PLSX-2510	10	25 cm	Fuchsia
PLSX-2510.5	10.5	25 cm	Light Fuchsia
PLSX-2511	11	25 cm	Yellow

*Data on file.



Understand. Innovate. Deliver.™

Merit Medical Systems, Inc.
1600 West Merit Parkway
South Jordan, Utah 84095
1.801.253.1600
1.800.35.MERIT

Merit Medical Europe, Middle East, & Africa (EMEA)
Amerikalaan 42, 6199 AE Maastricht-Airport
The Netherlands
+31 43 358 82 22

Merit Medical Ireland Ltd.
Parkmore Business Park West
Galway, Ireland
+353 (0) 91 703 733

Merit.com

Austria
0800 295 374

Belgium
0800 72 906 (Dutch)
0800 73 172 (Français)

Denmark
80 88 00 24

Finland
0800 770 586

France
0800 91 60 30

Germany
0800 182 0871

Ireland (Republic)
1800 553 163

Italy
800 897 005

Luxembourg
8002 25 22

Netherlands
0800 022 81 84

Norway
800 11629

Portugal
308 801 034

Russia
+7 495 221 89 02

Sweden
020 792 445

Switzerland
(Deutsch)
+41 225180252
(Français)
+41 225948000
(Italiano)
+41 225180035
UK
0800 973 115

MIN MEDICAL

MIN Medical GmbH | Franz Schalk Platz 9/2 | 1130 Wien | Austria
Mobil: +43 660 107 67 67 | Fax: +43 1 877 68 29 | Web: www.minmedical.at | Email: office@minmedical.at

401251001/B ID 041715