





CARDIAC INTERVENTION

PreludeSNAP™


Splittable Sheath Introducer



Easier to Break
requiring 66% less break
force than its predecessor



Seamless
sheath-to-dilator
transition



PreludeSNAP™

Splittable Sheath Introducer

The Prelude SNAP from Merit Medical features a **SEAMLESS** sheath-to-dilator transition to reduce vessel trauma and provide ease of insertion through tortuous anatomy. The Prelude SNAP is also **EASIER** to break, requiring 66% less break force than its predecessor, making it effortless to remove and helping ensure leads stay in place.




PRODUCT FEATURES

Robust Valve Design

Tight seal to protect against air embolism and blood loss.

Improved Grip

Expanded wings with non-slip, ridged surfaces facilitate a tighter, more assured grip.



Easier to Break

Requires 66% less break force than its predecessor.*

Seamless Transition

Smoother sheath-to-dilator transition to reduce vessel trauma and promote ease of insertion.

KIT CONTENTS

- Splittable Sheath
- Dilator
- Guide Wire 0.038" x 50 cm or 0.038" x 80 cm
- Merit Advance Needle 18G x 7cm
- 12 mL Syringe



Competitors' sheath-to-dilator transitions



Our seamless sheath-to-dilator transition

ORDERING INFORMATION

Catalog Number with Sideport	French Size	Sheath Length	Hub Color
PLS-1006	6	13 cm	Green
PLS-1007	7	13 cm	Orange
PLS-1008	8	13 cm	Blue
PLS-1008.5	8.5	13 cm	Light Blue
PLS-1009	9	13 cm	Black
PLS-1009.5	9.5	13 cm	Gray
PLS-1010	10	13 cm	Fuchsia
PLS-1010.5	10.5	13 cm	Light Fuchsia
PLS-1011	11	13 cm	Yellow
PLS-1012	12	13 cm	Brown
PLS-1012.5	12.5	13 cm	Light Brown
PLS-2506	6	25 cm	Green
PLS-2507	7	25 cm	Orange
PLS-2508	8	25 cm	Blue
PLS-2509	9	25 cm	Black
PLS-2510	10	25 cm	Fuchsia
PLS-2510.5	10.5	25 cm	Light Fuchsia
PLS-2511	11	25 cm	Yellow

Catalog Number without Sideport	French Size	Sheath Length	Hub Color
PLSX-1006	6	13 cm	Green
PLSX-1007	7	13 cm	Orange
PLSX-1008	8	13 cm	Blue
PLSX-1008.5	8.5	13 cm	Light Blue
PLSX-1009	9	13 cm	Black
PLSX-1009.5	9.5	13 cm	Gray
PLSX-1010	10	13 cm	Fuchsia
PLSX-1010.5	10.5	13 cm	Light Fuchsia
PLSX-1011	11	13 cm	Yellow
PLSX-1012	12	13 cm	Brown
PLSX-1012.5	12.5	13 cm	Light Brown
PLSX-2506	6	25 cm	Green
PLSX-2507	7	25 cm	Orange
PLSX-2508	8	25 cm	Blue
PLSX-2509	9	25 cm	Black
PLSX-2510	10	25 cm	Fuchsia
PLSX-2510.5	10.5	25 cm	Light Fuchsia
PLSX-2511	11	25 cm	Yellow

*Data on file.